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## ENGLISH

**1.** Open the pack and take out the NIO. Make sure that the NIO is free of all packaging parts.

**2.** Select one of the following injection sites:

**2A. Primary site for Intraosseous insertion: Proximal tibia**
Approximately 1 inch or 2 cm medially and 1/2 inch or 1 cm proximally to the tibial tuberosity.

**2B. Secondary site for Intraosseous insertion: Humeral Head**
Adduct the patient's hand and locate the greater tubercle next to the head of the humerus.

NOTE: To prevent accidental removal of the device, immobilize patient's arm after procedure is complete.

**3.** Disinfect the skin on the insertion site by following institutional protocols. Note that this is a two-handed procedure. Place your non-dormiant hand on the textured dots located on the lower part of the NIO and position the NIO at a 90 degrees angle to the skin at the insertion site. The non-dominant hand should maintain this position throughout this procedure.

**4.** Unlock the NIO by rotating the cap 90 degrees in either direction.

**5.** Place the palm of your dominant hand over the cap.

Press the device against the patient's skin and maintain downward pressure. While pressing down on the device, pull the trigger wings upwards. This action will activate the device.

**6.** Gently pull the NIO up in a rotary motion while holding the base of the needle stabilizer against the insertion site.

**7.** While holding the needle stabilizer and cannula in place, remove the stylet by pulling it up (use a twisting motion of the stylet if necessary). The keyhole notch on the distal end of the NIO can be used to assist in removing the stylet from the cannula. Place stylet into an appropriate biohazard container.

**8.** Connect a syringe and confirm secure fitting. If desired, aspirate bone marrow. Always confirm successful needle placement by flushing with up to 20 cc of fluid or per your protocol. It is recommended to use the NIO Fixation sticker to affix the NIO stabilizer. (A). Connect any standard system for infusion (B). Cover the insertion site with a sterile, occlusive dressing.

**9.** Removal instruction: Remove the cannula and needle stabilizer by twisting and pulling vertically. Dispose of the removed components in the appropriate biohazard container. Cover the insertion site with a sterile, occlusive, wound dressing.

**NIO DESIGN**

The NIO is a single use, semi-automatic, spring loaded IO device with dual safety features for enhanced caregiver and patient safety. After deployment, a unique stabilizer base firmly secures the needle in place.

**INDICATIONS FOR USE** NIO Adult (NIO-A)

The NIO is intended to provide intraosseous access in the proximal tibia, as an alternative to IV access during emergencies. The device is for use in adult patients only. The NIO is indicated for use in providing intraosseous access as an alternative to IV access during emergencies. Humeral head IO access is indicated when rapid fluid or pharmacological resuscitation is required and intravenous access is not possible. The device is for use in adult patients only.

**RECOMMENDATIONS** • Recommended penetration depth for Prox. Tibia / Humeral Head 2.5cm.

- Reassess IO site immediately after step 8 (confirm needle placement). It is recommended to frequently monitor the limb every 10 minutes for the first half hour or longer after beginning of drug administration.
- Use a pressure bag for optimal infusion rates.
- Continue to monitor extremity for complications on a regular basis, especially pre and post infusion. Before drug administration, the IO needle should be rechecked for placement and patency by confirming the needle is stable in the bone. Beside Doppler ultrasonography may be useful to ascertain IO placement and flow.
- For conscious patients, consider

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